Human Mitochondrial DNA Copy Number Assay Kit

(96 well assay suitable for 44 mitochondrial and 44 nuclear reactions, 4 positive controls and 4 blanks)

Catalog Number: MCN 1

Store at -20°C.

FOR RESEARCH USE ONLY

Introduction: This DNA analysis kit is for the determination of human mitochondrial DNA copy number, *in vivo* and *in vitro*, by the comparison of mitochondrial (mt) and nuclear (n) DNA measured by real-time PCR.

Kit Contents:

- 96 well PCR plate
- rtPCR reaction mix.
- Validated primers to quantify mitochondrial DNA (mtDNA).
- Validated primers to quantify nuclear DNA (nDNA)
- Positive control [1.825 ng/ μ l] (isolated total DNA from human MCF10A cells).

Not Included in Kit:

- DNA isolation Kit
- Nuclease-free water
- PCR Tubes and Caps

Thermal cycler program:

- Preprogram PCR machine for this profile:
 - a. 95°C, 10 min

(40 Cycles)

 - b. 95°C, 15 sec
 - c. 60°C, 60 sec

Real time PCR procedure: The following procedure is for each 20 μ L reaction. Increase all amounts proportionally according to the total number of tubes.

- Per PCR tube (20 μ L Rx), mix the following:
 - a. 1 μ L forward primer
 - b. 1 μ L reverse primer
 - c. 8 μ L sample contain genomic DNA/ 8 μ L of positive control
 - d. 10 μ L rtPCR reaction mix

Recommended concentration: Between 0.5 to 3.75ng/ μ L

Calculations: Mt copy number =


$$\Delta Ct1 = Ct (\text{mitochondria} - \text{control}) - Ct (\text{nucleus} - \text{control})$$

$\Delta Ct2 = Ct \text{ (mitochondria - experimental)} - Ct \text{ (nucleus - experimental)}$

$\Delta\Delta Ct = \text{Sample } \Delta Ct - \text{Average } \Delta Ct2 \text{ control.}$

mtDNA level change = $2^{-\Delta\Delta Ct}$

Total DNA isolated from human MCF10A cells

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
Plot of C_T versus DNA concentration

References

Horikawa, I. et al. (2017) $\Delta 133p53$ represses p53-inducible senescence genes and enhances the generation of human induced pluripotent stem cells.

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5442472/pdf/cdd201748a.pdf>

Guo, X. et al. (2020) STAMP: a multiplex sequencing method for simultaneous evaluation of mitochondrial DNA heteroplasmies and content.

<https://academic.oup.com/nargab/article/2/4/lqaa065/5930403>